




EUROArray HSV1/2 VZV



- Reliable direct detection of herpes simplex virus types 1 and 2 as well as varicella zoster virus in a single test
- Discrimination of pathogens even if clinical differentiation is difficult
- Suitable for swab samples from the skin, mouth or genital area

Technical data

Substrate	Single-stranded DNA probes, length: 15 to 50 nucleotides
Test procedure	DNA extraction / PCR (approx. 3 h) / hybridisation (60 min) / fully automated evaluation
Reagents	Ready for use
Controls	DNA-negative control and other integrated controls
Test kit format	20, 10 or 5 slides each containing 5 test fields, or 8 slides each containing 3 test fields
Order number	MN 2530-2005-1, MN 2530-1005-1, MN 2530-0505-1, MN 2530-0803-1

Clinical significance

Herpes simplex virus types 1 and 2 (HSV-1 and HSV-2) and varicella zoster virus (VZV) are human pathogenic viruses from the *Herpesviridae* family that occur worldwide.

HSV-1 and HSV-2, which are closely related, can cause different diseases, with either asymptomatic, mild or severe courses. The most common manifestations are herpes labialis (cold sore) and herpes genitalis, also rarely herpes encephalitis or neonatorum in severe cases. HSV-1 generally leads to cold sores. Genital infections can be caused by both virus types, but mainly by HSV-2. In primary infections the virus penetrates the mucosal cells, multiplies and is released in ulcers or inflammatory skin blisters. Infections persist lifelong, alternating between asymptomatic phases and acute infectious outbreaks. The virus is transmitted by blister fluid or genital secretions at skin or mucosal contact. HSV-1 is also transmitted via the saliva. Infection of the foetus in pregnancy is generally rare. However, if it occurs, it leads to intrauterine death, severe malformations or premature delivery in most cases.

The disease caused by primary infection with VZV is called chickenpox, a highly infectious childhood disease. VZV are transmitted by blister fluid, conjunctival fluid or saliva containing the virus. After first manifestation the virus establishes latency in sensory nerve cells and can subsequently reactivate to cause herpes zoster (shingles) as second manifestation. Typical symptoms are painful and/or itching skin exanthema, pains and sensory disorders. However, asymptomatic reactivation has also been described. In severe cases, complications such as pneumonia, hepatitis or post-zoster neuralgia can occur. Since reactivations are associated with a weakened cellular immune response, severe cases occur especially in elderly or immunosuppressed persons.

The EUROArray HSV1/2 VZV enables direct detection and differentiation between HSV-1, HSV-2 and VZV based on swab samples from the skin, mouth mucosa or genital mucosa of a patient. It also allows simple and quick discrimination of the pathogens, even if clinical differentiation is difficult. Moreover, it supports the differentiation of herpes infections from other dermatoses with similar symptoms.



Test principle


This test system is based on amplification of defined gene sections and subsequent detection via a hybridisation reaction with immobilised DNA probes in a microarray system. Isolated DNA from skin, oral or genital swabs of a patient is used as sample material. In the first reaction step, specific DNA sections from the pathogens present in the sample are amplified by polymerase chain reaction (PCR) using a multiplex primer system and, at the same time, labelled with a fluorescent dye. In a second step, the products are detected using an oligonucleotide microarray. The specific binding (hybridisation) of the fluorescence-labelled PCR product to the corresponding oligonucleotide probe is detected using a special microarray scanner. The EUROArrayScan software evaluates all spot signals automatically and deduces the test result.

DNA sample

PCR (polymerase
chain reaction)

DNA microarray
hybridisation

Fully automated
evaluation



Test performance

The PCRs are incubated in the thermocycler and then, using the TITERPLANE technique, incubated on EUROArray slides containing microarray BIOCHIPS. Scanning, evaluation and documentation are performed fully automatically using the EUROArray Scanner (incl. EUROArrayScan-Software).

Analytical sensitivity

The lower detection limit (LoD) of this test is 10 DNA copies/reaction for each of the three pathogens. The LOD is the minimum detection limit. Usually, fewer DNA copies of the pathogenic agent are detected.

Analytical specificity

The specificity of the test system is ensured by the primer and probe design, as well as the PCR and hybridisation conditions given in the test instructions. Using up to 2 million copies of genomic DNA of one of the pathogens to be detected, the test system did not show any cross reactivity with the other pathogens to be detected. Furthermore, no cross reactions were found for high copy numbers of nucleic acids of different other pathogens that occur in or on the human body.

Evaluation study

In an evaluation study, swab samples precharacterised by means of PCR-based reference tests were investigated with the EUROArray HSV1/2 VZV:

	HSV-1	HSV-2	VZV
Precharacterised positive samples,	110	101	114
thereof positive in EUROArray HSV1/2 VZV	106	101	107
Precharacterised negative samples,	191	200	136
thereof negative in EUROArray HSV1/2 VZV	190	199	133
Sensitivity	96.4%	100%	93.9%
Specificity	99.5%	99.5%	97.8%

Swabs from the skin, mouth and genital area as well as from the eye and anorectal region were investigated. The results from the analysis were 100% consistent with the precharacterisation in both sample groups.