



# A newly developed chemiluminescence immunoassay improves the detection of anti-PLA2R autoantibodies in primary membranous nephropathy

A. Kühnl<sup>1</sup>, C. Dährich<sup>1</sup>, S. Saschenbrecker<sup>1</sup>, I. Gunnarsson<sup>2,3</sup>, P. Ronco<sup>4,5,6</sup>,  
H. Debiec<sup>5,6</sup>, and W. Schlumberger<sup>1</sup>

<sup>1</sup>Institute for Experimental Immunology, affiliated to EUROIMMUN AG, Lübeck, Germany


<sup>2</sup>Division of Rheumatology, Department of Medicine, Karolinska Institutet, Stockholm, Sweden

<sup>3</sup>Rheumatology, Karolinska University Hospital, Stockholm, Sweden


<sup>4</sup>Department of Nephrology and Dialysis, Tenon Hospital, Paris, France

<sup>5</sup>Sorbonne Universités, UPMC University Paris 06, Paris, France


<sup>6</sup>Institut National de la Santé et de la Recherche Médicale (INSERM), Paris, France



Anti-PLA2R reactivity as determined by ChLIA, ELISA and RC-IFA in 309 sera (155 pMN, 154 disease controls). Percentages indicate the agreement between assays.



Correlation between ChLIA and ELISA for the detection of anti-PLA2R autoantibodies in sera from 155 pMN patients. Dashed lines represent cut-off values for positivity.



Correlation between ChLIA and RC-IFA for the detection of anti-PLA2R autoantibodies in sera from 155 pMN patients. Dashed lines represent cut-off values for positivity.

## Introduction

Autoantibodies against the **M-type phospholipase A2 receptor (PLA2R)** are important markers in the diagnosis and monitoring of **primary membranous nephropathy (pMN)**. For the detection of anti-PLA2R autoantibodies, a standardised recombinant **cell-based indirect immunofluorescence assay (RC-IFA)** and ELISA are widely used, the former providing higher sensitivity but lacking a finely graduated quantification of

antibody titers. In this study, we evaluated the diagnostic performance characteristics of a novel standardised **chemiluminescence immunoassay (ChLIA)** by comparison with the established anti-PLA2R test systems.

## Methods

Sera from 155 biopsy-proven pMN patients and 154 disease controls were analysed for autoantibodies against PLA2R by the novel ChLIA as well as by ELISA and RC-IFA.

## Results

The clinical sensitivity of the ChLIA (83.9%) was higher compared to ELISA (73.5%) and equalled that of RC-IFA (83.2%), at similar specificities ( $\geq 99.4\%$ ). Among ELISA-negative pMN samples, ChLIA and RC-IFA yielded positive results in 39.0% and 36.6% of cases, respectively. The overall qualitative agreement amounted to 94.5% (ChLIA vs. ELISA) and 99.4% (ChLIA vs. RC-IFA).

## Conclusion

The novel anti-PLA2R ChLIA outperforms the ELISA in detecting pMN patients and demonstrates almost perfect agreement with RC-IFA. It thus presents a promising alternative tool for accurate anti-PLA2R testing in routine diagnostic settings, with the advantage of rapid turnaround times and fully automated random-access processing.

ROC curve analysis	Anti-PLA2R positive		
	ChLIA	ELISA	RC-IFA
AUC	0.899	0.927	0.916
Max. sum of sensitivity and specificity (cut-off)	183.9% (9.1)	181.3% (3.0)	183.2% (1.5)
Sensitivity at 98% specificity (cut-off)	85.2% (7.8)	81.9% (5.8)	NA
Sensitivity at 99% specificity (cut-off)	84.5% (9.1)	80.7% (7.3)	NA
Sensitivity at 100% specificity (cut-off)	83.2% (15.5)	76.1% (15.2)	83.2% (1.5)

Panel	n	Anti-PLA2R positive		
		ChLIA	ELISA	RC-IFA
Primary MN	155	130	114	129
Sensitivity	155	83.9%	73.5%	83.2%
Secondary MN	6	0	0	0
IgA nephropathy	10	0	0	0
FS glomerular sclerosis	10	0	0	0
MP glomerulonephritis	10	0	0	0
Minimal change disease	17	1	0	0
Lupus nephritis class I-V	33	0	0	0
Lupus nephritis class V	34	0	0	0
System. lupus erythem.	34	0	0	0
Specificity	154	99.4%	100%	100%